


GOAT ANTI-BIRC3 (AA 109 -118) ANTIBODY

SKU: EB10370



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	BIRC3 TNFR2-TRAF signaling complex protein RNF49 MIHC MALT2 HIAP1 HAIP1 c-IAP2 CIAP2 baculoviral IAP repeat-containing protein 3 baculoviral IAP repeat-containing 3 apoptosis inhibitor 2 API2 AIP1
Accession ID	NP_001156.1
Blocking Peptide	EBP10370
Immunogen	Peptide with sequence C-SQPTFPSSVT, from the internal region of the protein sequence according to NP_001156.1.
Product Comments	Reported variants represent identical protein: NP_892007.1, NP_001156.1
Peptide Sequence	C-SQPTFPSSVT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	330
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Tonsil. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 70kDa band observed in lysates of cell line HeLa (calculated MW of 68.4kDa according to NP_001156.1). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

