


GOAT ANTI-BLM ANTIBODY

SKU: EB09998



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	BLM RECQL3 RECQL2 RECQ2 MGC131620 MGC131618 MGC126616 BS Bloom syndrome, RecQ helicase-like Bloom syndrome protein
Names	
Accession ID	NP_000048.1
Blocking Peptide	EBP09998
Immunogen	Peptide with sequence C-ESVVRESKSFTRST, from the C Terminus of the protein sequence according to NP_000048.1.
Product Comments	Wrong peptide for different protein
Peptide Sequence	C-ESVVRESKSFTRST
Shipping Instructions	Refrigerated
Predicted Species	Human, Rat
Reactive Species	Human
Human Gene ID	641
Rat Gene ID	113959
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 160kDa band observed in nuclear lysates of cell lines HeLa and Jurkat (calculated MW of 159kDa according to NP_000048.1). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

