



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600


GOAT ANTI-BMP4 ANTIBODY

SKU: EB12544



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600







Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

bone morphogenetic protein 2B|BMP-4|BMP-2B|ZYME|OFC11|MCOPS6|BMP2B1|BMP2B|bone morphogenetic

protein 4|BMP4

Accession

Names

NP_001193.2

Blocking

Peptide

EBP12544

Immunogen

Peptide with sequence C-SPKHHSQRARK, from the internal region of the protein sequence according to

NP 001193.2.

Product

The immunizing peptide represents the N terminus of the mature protein. Reported variants represent

Comments

identical protein: NP_570911.2, NP_001193.2, NP_570912.2.

Peptide

Sequence

C-SPKHHSQRARK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

Reactive

Species Human

652

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:64000.

NP_001193.2). Recommended concentration: 0.3-1µg/ml.

Limit

Western Blot

Approx 40-45kDa band observed in lysates of cell lines A431 and HeLa (calculated MW of 46.6kDa according to

Application

Type

Pep-ELISA, WB

GALLERY IMAGES





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

