

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-BNC ANTIBODY

SKU: EB05052

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias zinc finger protein basonuclin|basonuclin|BNC

Names

Accession NP 001708 ID

Blocking

EBP05052 Peptide

Immunogen NP_001708. Peptide with sequence C-NLHKSLASSPSHLQ, from the C Terminus of the protein sequence according to

Peptide

C-NLHKSLASSPSHLQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse **Species**

Human

646 **Gene ID**

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Western

Western blot: No signal obtained yet but low background observed in Jurkat and HeLa lysates at upto 2µg/ml. Blot