

GOAT ANTI-BNC ANTIBODY

SKU: EB05052



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	zinc finger protein basonuclin basonuclin BNC
Accession ID	NP_001708
Blocking Peptide	EBP05052
Immunogen	Peptide with sequence C-NLHKSLASSPSHLQ, from the C Terminus of the protein sequence according to NP_001708.
Peptide Sequence	C-NLHKSLASSPSHLQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	646
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Jurkat and HeLa lysates at upto 2µg/ml.