

GOAT ANTI-BRADYKININ RECEPTOR B1 ANTIBODY

SKU: EB08800



250kDa
150kDa
100kDa
75kDa

50kDa
37kDa

25kDa
20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	bradykinin receptor 1 bradykinin B1 receptor BK-1 receptor BRADYB1 BKR1 BKB1R B1R B1BKR BDKRB1 bradykinin receptor B1
Accession ID	NP_000701.2
Blocking Peptide	EBP08800
Immunogen	Peptide with sequence KVVWELYKQCTPK, from the internal region of the protein sequence according to NP_000701.2.
Peptide Sequence	KVVWELYKQCTPK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	623
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx 35kDa band observed in lysates of cell lines Daudi, K562 and U937 calculated MW of 40.5kDa according to NP_000701.2). This molecular weight is routinely observed by other sources. Recommended concentration: 2-4µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

