

GOAT ANTI-BRDG1 ANTIBODY

SKU: EB05108



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	STAP-1 BCR downstream signaling 1 BRDG1
Names	
Usage Summary	Please note that we observe approximately 40kDa band instead of the predicted 32kDa. Ohya et al., (PNAS, Oct 1999; 96: 11976 - 11981) also observe a band around 40kDa using rabbit anti-BRDG1.
Blocking Peptide	EBP05108
Immunogen	Peptide with sequence C-GRSEKLKKNPHIA, from the C Terminus of the protein sequence according to .
Peptide Sequence	C-GRSEKLKKNPHIA
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	26228
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 40kDa bands seen in Raji lysate. Recommended for use at 1-3µg/ml.