

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-BRDG1 ANTIBODY

SKU: EB05108

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias STAP-1|BCR downstream signaling 1|BRDG1

Names

Usage Please note that we observe approximately 40kDa band instead of the predicted 32kDa. Ohya et al., (PNAS,

Summary Oct 1999; 96: 11976 - 11981) also observe a band around 40kDa using rabbit anti-BRDG1.

Blocking

Peptide

Immunogen Peptide with sequence C-GRSEKLKKNPHIA, from the C Terminus of the protein sequence according to .

Peptide

C-GRSEKLKKNPHIA Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human

Species

Human

Gene ID

26228

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Western

Western blot: Approx 40kDa bands seen in Raji lysate. Recommended for use at 1-3µg/ml. **Blot**