

GOAT ANTI-BRUNO-LIKE 5 ANTIBODY

SKU: EB08045



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	bruno-like 5, RNA binding protein RNA-binding protein BRUNOL-5 CUG-BP and ETR-3 like factor 5 Bruno (Drosophila) -like 5, RNA binding protein CELF5 BRUNOL-5 bruno-like 5, RNA binding protein (Drosophila) BRUNOL5
Accession ID	NP_068757.2
Blocking Peptide	EBP08045
Immunogen	Peptide with sequence C-EPPGGQPDGMKD, from the internal region (near the N Terminus) of the protein sequence according to NP_068757.2.
Peptide Sequence	C-EPPGGQPDGMKD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	60680
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx 55kDa band observed in Human Heart lysates (calculated MW of 52.4kDa according to NP_068757.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

