


GOAT ANTI-C16ORF57 (AA154-167) ANTIBODY

SKU: EB11248



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	C16orf57 PN hypothetical protein LOC79650 HVSL1 HVSL motif containing 1 FLJ13154 chromosome 16 open reading frame 57
Accession ID	NP_078874.2; NP_001182231.1
Blocking Peptide	EBP11248
Immunogen	Peptide with sequence C-TANQVKIYTNQEK, from the internal region of the protein sequence according to NP_078874.2; NP_001182231.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_078874.2; NP_001182231.1).
Peptide Sequence	C-TANQVKIYTNQEK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Pig
Reactive Species	Human, Mouse, Rat
Human Gene ID	79650
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 35+28kDa bands observed in Human, Mouse and Rat Skin lysates (calculated MW of 30.3kDa according to NP_078874.2 and 28.1kDa according to NP_001182231.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

