

GOAT ANTI-C20ORF75 ANTIBODY

SKU: EB08262



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias / Names	LRRN4 neuronal leucine-rich repeat protein 4 hypothetical protein LOC164312 dJ1056H1.1 NLRR4 MGC25027 FLJ23994 RP5-967N21.9 chromosome 20 open reading frame 75 C20orf75
Accession ID	NP_689824.1
Blocking Peptide	EBP08262
Immunogen	Peptide with sequence C-YKNPAFDDYP, from the C Terminus of the protein sequence according to NP_689824.1.
Peptide Sequence	C-YKNPAFDDYP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	164312
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA