

GOAT ANTI-CASPASE 12 (AA232-244) ANTIBODY

SKU: EB11672

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	CASP12 CASP-12 CASP12P1 caspase 12 (gene/pseudogene) caspase 12 pseudogene 1 inactive caspase-12 OTTHUMP00000207026 OTTHUMP00000207028 OTTHUMP00000207029 OTTHUMP00000207031 OTTHUMP00000207032 OTTHUMP00000219058
Accession ID	NP_001177945.1
Blocking Peptide	EBP11672
Immunogen	Peptide with sequence C-TDSGKAGADTHGR, from the internal region of the protein sequence according to NP_001177945.1.
Peptide Sequence	C-TDSGKAGADTHGR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	100506742
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 35kDa band observed in Human Heart lysates (calculated MW of 38.9kDa according to NP_001177945.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

