

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-CASPASE 3 ANTIBODY**

**SKU:** EB05257

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

SREBP cleavage activity 1|cysteine protease CPP32|PARP cleavage **Alias** 

**Names** 

**Accession** NP\_004337.2; NP\_116786.1 ID

Peptide with sequence RGTELDCGIETD, from the internal region of the protein sequence according to **Immunogen** 

NP 004337.2; NP 116786.1.

**Peptide** 

**RGTELDCGIETD** Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

protease|apopain|Yama|APOPAIN|CPP32B|SCA-1|CPP32|caspase 3, apoptosis-related cysteine protease|CASP3

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Hamster, Dog, Sheep, Pig **Species** 

Reactive

Human **Species** 

Human

836

**Gene ID** 

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png Grade

Western

Western blot: Approx 40kDa band seen in A431 cell lysates. Recommended for use at 1-3µg/ml. **Blot**