

GOAT ANTI-CASPASE 3 ANTIBODY

SKU: EB05257



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SREBP cleavage activity 1 cysteine protease CPP32 PARP cleavage protease apopain Yama APOPAIN CPP32B SCA-1 CPP32 caspase 3, apoptosis-related cysteine protease CASP3
Accession ID	NP_004337.2; NP_116786.1
Immunogen	Peptide with sequence RGTELDGIEDT, from the internal region of the protein sequence according to NP_004337.2; NP_116786.1.
Peptide Sequence	RGTELDGIEDT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Hamster, Dog, Sheep, Pig
Reactive Species	Human
Human Gene ID	836
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 40kDa band seen in A431 cell lysates. Recommended for use at 1-3µg/ml.