

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CASPASE 6 ANTIBODY

SKU: EB12166



Email: customerservice@vectorlabs.com

Telephone: <u>(650)</u> 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

CASP6|caspase 6, apoptosis-related cysteine peptidase|MCH2|apoptotic protease MCH-2|caspase 6, apoptosis-

related cysteine protease|caspase-6

Accession

Names

NP_001217.2; NP_116787.1

Blocking

EBP12166 **Peptide**

Peptide with sequence C-QTEKLDTNITEVD, from the internal region of the protein sequence according to **Immunogen**

NP_001217.2; NP_116787.1.

Product Comments

This antibody is expected to recognize both reported isoforms (NP 001217.2; NP 116787.1).

Peptide

Sequence

C-QTEKLDTNITEVD

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Reactive

Human

Species Human

839 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

ELISA

Detection

Antibody detection limit dilution 1:128000.

Limit

Western Approx 33kDa band observed in Human Colon lysates (calculated MW of 33.3kDa according to NP_001217.2).

Recommended concentration: 0.3-1µg/ml. Blot

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa