

GOAT ANTI-CASPASE 6 ANTIBODY

SKU: EB12166

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	CASP6 caspase 6, apoptosis-related cysteine peptidase MCH2 apoptotic protease MCH-2 caspase 6, apoptosis-related cysteine protease caspase-6
Accession ID	NP_001217.2; NP_116787.1
Blocking Peptide	EBP12166
Immunogen	Peptide with sequence C-QTEKLDNITEVD, from the internal region of the protein sequence according to NP_001217.2; NP_116787.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_001217.2; NP_116787.1).
Peptide Sequence	C-QTEKLDNITEVD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	839
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 33kDa band observed in Human Colon lysates (calculated MW of 33.3kDa according to NP_001217.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

