

GOAT ANTI-CBL ANTIBODY

SKU: EB05288



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa

A Western blot image showing a single, prominent dark band at approximately 150kDa. The molecular weight markers are indicated on the left side of the blot, ranging from 250kDa at the top to 25kDa at the bottom. The band is located between the 150kDa and 100kDa markers.

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	oncogene CBL2 Cas-Br-M ecotropic retroviral transforming sequence (Oncogene CBL2) Cas-Br-M (murine) ecotropic retroviral transforming sequence RNF55 C-CBL CBL2 CBL
Accession ID	NP_005179.2
Blocking Peptide	EBP05288
Immunogen	Peptide with sequence C-REFVSISSPAHVAT, from the C Terminus of the protein sequence according to NP_005179.2.
Peptide Sequence	C-REFVSISSPAHVAT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Mouse
Human Gene ID	867
Mouse Gene ID	12402
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 120kDa band seen in 3T3 lysate. Recommended for use at 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

