

GOAT ANTI-CCL21 ANTIBODY

SKU: EB10008

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	CCL21 TCA4 small inducible cytokine A21 SLC secondary lymphoid tissue chemokine SCYA21 MGC34555 exodus-2 Efficient Chemoattractant for Lymphocytes ECL CKb9 chemokine (C-C motif) ligand 21 beta chemokine exodus-2 6Ckine
Accession ID	NP_002980.1
Blocking Peptide	EBP10008
Immunogen	Peptide with sequence QGCRKDRGASKTG, from the internal region of the protein sequence according to NP_002980.1.
Product Comments	This antibody may cross-react with mouse CCL21 (GeneID: 18829) and rat CCL21 (GeneID: 298006).
Peptide Sequence	QGCRKDRGASKTG
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	6366
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx 16kDa band observed in Human Lymph Node lysates (calculated MW of 14.6kDa according to NP_002980.1). Recommended concentration: 0.03-0.1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

