

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CCL3L1 ANTIBODY

SKU: EB09000

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / small inducible cytokine A3-like

1|SCYA3L1|SCYA3L|MGC12815|MGC104178|LD78BETA|LD78|G0S19-2|D17S1718|464.2|chemokine (C-C motif) Alias

Names ligand 3-like 1|CCL3L1

Accession

NP 066286.1 ID

Blocking

EBP09000 **Peptide**

Peptide with sequence CNQVLSAPLAADT, from the internal region of the protein sequence according to Immunogen

NP 066286.1.

Peptide

CNQVLSAPLAADT Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

6349 Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection

Antibody detection limit dilution 1:8000.

Limit

Western Blot

CCL3L1 expression requires activation.

Application

Pep-ELISA **Type**