

GOAT ANTI-CCL3L1 ANTIBODY

SKU: EB09000



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	small inducible cytokine A3-like
Names	1 SCYA3L1 SCYA3L MGC12815 MGC104178 LD78BETA LD78 G0S19-2 D17S1718 464.2 chemokine (C-C motif) ligand 3-like 1 CCL3L1
Accession ID	NP_066286.1
Blocking Peptide	EBP09000
Immunogen	Peptide with sequence CNQVLSAPLAADT , from the internal region of the protein sequence according to NP_066286.1.
Peptide Sequence	CNQVLSAPLAADT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	6349
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	CCL3L1 expression requires activation.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)