


## GOAT ANTI-CCP5 (MOUSE) ANTIBODY

SKU: EB11040



## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Synonym / Alias Names</b> | Agbl5  ATP/GTP binding protein-like 5  4930455N08  9430057O19Rik  CCP5  ATP/GTP-binding protein-like 5  OTTMUSP00000025344  OTTMUSP00000025346  OTTMUSP00000025347  OTTMUSP00000025350  OTTMUSP00000025351  carboxypeptidase 5, cytosolic  4930455N08  9430057O19Rik  ATP/GTP binding protein-like 5  ATP/GTP-binding protein-like 5  carboxypeptidase 5, cytosolic  CCP5  cytosolic carboxypeptidase-like protein 5  OTTMUSP00000025344  OTTMUSP00000025346  OTTMUSP00000025347  OTTMUSP00000025348  OTTMUSP00000025350  OTTMUSP00000025351 Agbl5 |
| <b>Usage Summary</b>         | <strong>Flow Cytometry:</strong> Flow cytometric analysis of NIH3T3 cells. Recommended concentration: 10ug/ml.                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Accession ID</b>          | NP_001041657.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Blocking Peptide</b>      | EBP11040                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Immunogen</b>             | Peptide with sequence CFSKPEEAGSHVE, from the internal region of the protein sequence according to NP_001041657.1.                                                                                                                                                                                                                                                                                                                                                                                                                                 |

|                              |                                                                                                                                                                                                                                   |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Comments</b>      | This antibody is expected to recognize reported isoform 3 (NP_001041657.1) only.                                                                                                                                                  |
| <b>Peptide Sequence</b>      | CFSKPEEAGSHVE                                                                                                                                                                                                                     |
| <b>Shipping Instructions</b> | Refrigerated                                                                                                                                                                                                                      |
| <b>Predicted Species</b>     | Mouse                                                                                                                                                                                                                             |
| <b>Reactive Species</b>      | Mouse                                                                                                                                                                                                                             |
| <b>Mouse Gene ID</b>         | 231093                                                                                                                                                                                                                            |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |
| <b>ELISA Detection Limit</b> | Antibody detection limit dilution 1:16000.                                                                                                                                                                                        |
| <b>Application Type</b>      | Pep-ELISA, FC                                                                                                                                                                                                                     |

## GALLERY IMAGES

