

GOAT ANTI-CCRN4L / NOCTURNIN ANTIBODY

SKU: EB08119



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	nocturnin CCR4-like (carbon catabolite repression 4, <i>S.cerevisiae</i>) CCR4 carbon catabolite repression 4-like NOC MGC78549 MGC4120817 MGC142060 MGC142054 CCR4L CCR4 carbon catabolite repression 4-like (<i>S. cerevisiae</i>) CCRN4L
Accession ID	NP_036250.2
Immunogen	Peptide with sequence HLSLVCDFSFTESD, from the C Terminus of the protein sequence according to NP_036250.2.
Peptide Sequence	HLSLVCDFSFTESD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	25819
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Western Blot: Approx 37kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 48.2kDa according to NP_036250.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Garbarino-Pico et al, jRNA. 2007 May;13(5):745-55; PMID: 17400819). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

