Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CCRN4L / NOCTURNIN ANTIBODY

SKU: EB08119

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Alias

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

nocturnin|CCR4-like (carbon catabolite repression 4, S.cerevisiae)|CCR4 carbon catabolite repression 4-Synonym /

like|NOC|MGC78549|MGC4120817|MGC142060|MGC142054|CCR4L|CCR4 carbon catabolite repression 4-like

Names (S. cerevisiae)|CCRN4L

Accession

NP 036250.2 ID

Peptide with sequence HLSLVCDFSFTEESD, from the C Terminus of the protein sequence according to **Immunogen**

NP 036250.2.

Peptide

HLSLVCDFSFTEESD Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human

Species Human

25819

Gene ID

Product

Grade

ELISA

Limit

Detection Antibody detection limit dilution 1:4000.

Western

Blot

Western Blot: Approx 37kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 48.2kDa according to NP 036250.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Garbarino-Pico et al, jRNA. 2007 May;13(5):745-55; PMID: 17400819). Recommended

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB **Type**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600