

GOAT ANTI-CD20 / MS4A1 (C TERMINUS), BIOTINYLATED ANTIBODY

SKU: EB05140-B

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MGC3969 leukocyte surface antigen Leu-16 CD20 receptor CD20 antigen B-lymphocyte cell-surface antigen B1 S7 MS4A2 LEU-16 CVID5 CD20 Bp35 B1 membrane-spanning 4-domains, subfamily A, member 1 MS4A1
Accession ID	NP_068769.2
Blocking Peptide	EBP05140-B
Immunogen	Peptide with sequence CQDQESSPIENDSSP., from the C Terminus of the protein sequence according to NP_068769.2.
Product Comments	Reported variants represent identical protein: NP_068769.2, NP_690605.1
Peptide Sequence	CQDQESSPIENDSSP.
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	931
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 35kDa band observed in lysates of cell line Daudi (calculated MW of 33.1kDa according to NP_068769.2). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa