

GOAT ANTI-CD25 / IL2RA (AA117-129) ANTIBODY

SKU: EB11856



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	TCGFR TAC antigen p55 OTTHUMP00000019032 OTTHUMP00000019031 OTTHUMP00000019030 interleukin-2 receptor subunit alpha interleukin 2 receptor, alpha IL-2-RA IL2-RA IL2RA IL-2R subunit alpha IL2R IL-2 receptor subunit alpha IDDM10 CD25
Accession ID	NP_000408.1
Blocking Peptide	EBP11856
Immunogen	Peptide with sequence DQASLPGHCREPP, from the internal region of the protein sequence according to NP_000408.1.
Peptide Sequence	DQASLPGHCREPP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	3559
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments in Human Colorectal cancer, Cerebellum, Heart and CACO-2 lysates gave no specific signal but low background (at antibody concentration up to 2µg/ml).
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)