


GOAT ANTI-CD81 ANTIBODY

SKU: EB07017



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	target of antiproliferative antibody 1 CD81 antigen 26 kDa cell surface protein TAPA-1 TAPA1 S5.7 HGNC:1701 CD81 antigen (target of antiproliferative antibody 1) CD81
Usage Summary	Flow Cytometry: Flow cytometric analysis of Jurkat cells. Recommended concentration: 10ug/ml.
Accession ID	NP_004347.1; NP_001284578.1
Blocking Peptide	EBP07017
Immunogen	Peptide with sequence NLFKEDCHQKIDD, from the internal region of the protein sequence according to NP_004347.1; NP_001284578.1.
Peptide Sequence	NLFKEDCHQKIDD

Shipping Instructions	Refrigerated
Predicted Species	Human, Rat, Dog
Reactive Species	Human
Human Gene ID	975
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, FC

GALLERY IMAGES

