





Telephone: (650) 697-3600

GOAT ANTI-CDO ANTIBODY

SKU: EB09472

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

surface glycoprotein, Ig superfamily member|cell adhesion molecule-related/down-regulated by **Alias** oncogenes|ORCAM|MGC111524|CDO|Cdon homolog (mouse)|CDON

Names

Accession NP 058648.4 ID

Blocking

EBP09472 **Peptide**

Peptide with sequence C-KNRQQNTIQKYDP, from the internal region of the protein sequence according to Immunogen

NP 058648.4.

Peptide C-KNRQQNTIQKYDP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Dog, Cow **Species**

Human

50937 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA **Type**