

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-CDT1 / DUP ANTIBODY**

**SKU:** EB05607



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

DNA replication factor, Double parked, Drosophila, homolog of Double parked RIS2 chromatin licensing and DNA replication factor 1|DUP|CDT1

Names

**Accession** 

NP\_112190.2 ID

**Blocking** 

EBP05607 **Peptide** 

Peptide with sequence C-ARLAHQTRAEEGL, from the C Terminus of the protein sequence according to **Immunogen** 

NP 112190.2.

**Peptide** Sequence

C-ARLAHOTRAEEGL

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive

Human **Species** 

Human

81620 Gene ID

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:64000.

Limit

Western Approx 60kDa band observed in Human Ovary lysates (calculated MW of 60.4kDa according to NP\_112190.2).

**Blot** Recommended concentration: 1-3µg/ml.

**Application** 

Pep-ELISA, WB Type

## **GALLERY IMAGES**





Telephone: (650) 697-3600

250kDa 150kDa 100kDa

ederest

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa