

GOAT ANTI-CDY2 ANTIBODY

SKU: EB05778



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	testis-specific chromodomain protein Y protein 2 Y chromosome chromodomain protein 2A CDY2 CDY chromodomain protein, Y-linked, 2A chromodomain protein, Y chromosome, 2 CDY2A
Accession ID	NP_004816.1
Blocking Peptide	EBP05778
Immunogen	Peptide with sequence C-SMLKYVENKIDEF, from the C Terminus of the protein sequence according to NP_004816.1.
Peptide Sequence	C-SMLKYVENKIDEF
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	9426
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments gave no signal but low background in Hela and 293 whole-cell lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA