

GOAT ANTI-CEBPB (AA68-81) ANTIBODY

SKU: EB12113

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	CEBPB CCAAT/enhancer binding protein (C/EBP), beta C/EBP-beta CRP2 IL6DBP LAP NF-IL6 TCF5 CCAAT/enhancer-binding protein beta TCF-5 interleukin 6-dependent DNA-binding protein liver-enriched transcriptional activator protein nuclear factor NF-IL6 nuclear factor of interleukin 6 transcription factor 5
Accession ID	NP_005185.2
Blocking Peptide	EBP12113
Immunogen	Peptide with sequence C-DHERAIDFSPYLEP, from the internal region of the protein sequence according to NP_005185.2.
Peptide Sequence	C-DHERAIDFSPYLEP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Reactive Species	Human, Mouse, Pig
Human Gene ID	1051
Mouse Gene ID	12608
Rat Gene ID	24253
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 35kDa band observed in Human Skin and Pig Lung lysates, and in lysates of cell line NIH3T3 (calculated MW of 36.1kDa according to Human NP_005185.2, 31.5kDa according to Mouse NP_034013.1 and 36.2kDa according to Pig NP_001186818.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

