

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CEBPB (AA68-81) ANTIBODY

SKU: EB12113



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

CEBPB|CCAAT/enhancer binding protein (C/EBP), beta|C/EBP-beta|CRP2|IL6DBP|LAP|NF-IL6|TCF5|CCAAT/enhancer-binding protein beta|TCF-5|interleukin 6-dependent DNA-binding protein|liver-

enriched transcriptional activator protein|nuclear factor NF-IL6|nuclear factor of interleukin 6|transcription **Names** factor 5

Accession

NP 005185.2

Blocking

EBP12113 **Peptide**

Peptide with sequence C-DHERAIDFSPYLEP, from the internal region of the protein sequence according to **Immunogen**

NP 005185.2.

Peptide

C-DHERAIDFSPYLEP **Sequence**

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Pig, Cow

Reactive

Human, Mouse, Pig **Species**

Human

1051 Gene ID

Mouse **Gene ID**

12608

Rat Gene ID 24253

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

FLISA

Detection Antibody detection limit dilution 1:128000.

Limit

Blot

Approx 35kDa band observed in Human Skin and Pig Lung lysates, and in lysates of cell line NIH3T3

Western (calculated MW of 36.1kDa according to Human NP_005185.2, 31.5kDa according to Mouse NP_034013.1 and

36.2kDa according to Pig NP_001186818.1). Recommended concentration: 1-3µg/ml.

Application

Pep-ELISA, WB **Type**





Telephone: (650) 697-3600

GAL	LERY	IMAGES

LIVI IIIAGLS			
250kDa 150kDa	250kDa	Α	В
100kDa	150kDa 100kDa		
75kDa	75kDa		
	50kDa		
50kDa	37kDa		_
37kDa			
	25kDa		
	20kDa		

25kDa

20kDa

15kDa

15kDa