


GOAT ANTI-CENTA2 ANTIBODY

SKU: EB06121



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names cent-b, centaurin beta|ArfGAP with dual PH domains 2|centaurin-alpha 2 protein|HSA272195|centaurin, alpha 2|CENTA2|ADAP2

Accession ID NP_060874.1

Blocking Peptide EBP06121

Immunogen Peptide with sequence C-RLTASTESGRSSR, from the C Terminus of the protein sequence according to NP_060874.1.

Peptide Sequence C-RLTASTESGRSSR

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human
Reactive Species	Human
Human Gene ID	55803
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Placenta shows strong staining of cytoplasm in syncytiotrophoblasts. Recommended concentration, 5-10µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

