

GOAT ANTI-CENTG2 / AGAP1 ANTIBODY

SKU: EB06030



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	MGC71657 ArfGAP with GTPase domain, ankyrin repeat and PH domain 1 Arf GAP with GTP-binding protein-like, ANK repeat and PH domains 1 GTP-binding and GTPase-activating protein 1 centaurin, gamma
Names	2 KIAA1099 GGAP1 AGAP1 CENTG2
Accession ID	NP_001032208.1; NP_055729.2
Blocking Peptide	EBP06030
Immunogen	Peptide with sequence C-NRNNSSGRVPTII, from the C Terminus of the protein sequence according to NP_001032208.1; NP_055729.2.
Product Comments	This antibody is expected to recognise both reported isoforms.
Peptide Sequence	C-NRNNSSGRVPTII
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog
Human Gene ID	116987
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Preliminary experiments gave no signal but low background in human kidney and brain lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

