Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CENTG2 / AGAP1 ANTIBODY

SKU: EB06030

SPECIFICATIONS

Unit Size 100 µg

Synonym / MGC71657|ArfGAP with GTPase domain, ankyrin repeat and PH domain 1|Arf GAP with GTP-binding protein-like,

Alias ANK repeat and PH domains 1|GTP-binding and GTPase-activating protein 1|centaurin, gamma

Names 2|KIAA1099|GGAP1|AGAP1|CENTG2

Accession

ID **Blocking**

EBP06030 **Peptide**

Peptide with sequence C-NRNNSSGRVPTII, from the C Terminus of the protein sequence according to **Immunogen**

NP 001032208.1; NP 055729.2.

NP_001032208.1; NP_055729.2

Product

This antibody is expected to recognise both reported isoforms. Comments

Peptide

C-NRNNSSGRVPTII Sequence

Shipping

Instructions Refrigerated

Predicted

Human, Dog **Species**

Human 116987

Gene ID

Product

Grade

FI ISA

Detection Antibody detection limit dilution 1:4000.

Limit

Blot

Western

Preliminary experiments gave no signal but low background in human kidney and brain lysates at up to 1µg/ml.

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

We would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?

Application

Pep-ELISA **Type**