


GOAT ANTI-CENTROMERE PROTEIN F / MITOSIN ANTIBODY

SKU: EB05622



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names mitosin|centromere protein F|hcp-1|RP11-262H5.1|centromere protein F, 350/400ka (mitosin)|cell-cycle-dependent 350K nuclear protein|centromere protein F (400kD)|CENP-F kinetochore protein|AH antigen|mitosin|PRO1779|CENF|centromere protein F (350/400kD, mitosin)|CENPF

Accession ID AAA82889.1; AAA82935.1; AAA86889.1

Blocking Peptide EBP05622

Immunogen Peptide with sequence GLESKGSENCKVQ, from the C Terminus of the protein sequence according to AAA82889.1; AAA82935.1; AAA86889.1.

Product Comments The immunizing peptide shows a K to N mismatch with the corresponding sequence of the NCBI-curated isoform NP_057427.3.

Peptide Sequence GLESKGSENCKVQ

Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Reactive Species	Human
Human Gene ID	1063
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Testis. Recommended concentration: 8µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

