


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CEP290 ANTIBODY

SKU: EB09122



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias

prostate cancer antigen T21|nephrocytsin-6|monoclonal antibody 3H11 antigen|CTCL tumor antigen se2-2|3H11Ag, BBS14, FLJ13615, FLJ21979, JBTS5, JBTS6, KIAA0373, LCA10, MKS4, NPHP6, SLSN6,

Names

rd16|centrosomal protein 290kDa|CEP290

Usage

Immunocytochemistry: Formaldehyde-fixed cell lines derived from Prostate Cancer and

Summary

from Breast Cancer show typical centrosomal staining. Recommended concentration: 2-3ug/ml.

Accession ID

NP 079390.3

EBP09122

Blocking

Peptide

Immunogen

Peptide with sequence C-RNSKHLKQQQYRAEN, from the internal region of the protein sequence according to NP 079390.3.

Product

Despite similarity of the immunizing peptide with the mouse protein, this antibody was demonstrated to weak

Comments for mouse Cep290.

Peptide Sequence

C-RNSKHLKQQQYRAEN

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human

Species Reactive

Human **Species**





Telephone: (650) 697-3600



Human Gene ID

80184

Mouse

216274

Gene ID Product

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite_medium.png$

IHC Results Paraffin embedded Human Prostate. Recommended concentration: 5-7µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:128000.

Limit

Application Pep-ELISA, IHC, ICC

Type

SELECTED REFERENCES

[{"pmid": 0, "intro": "", "title": "", "author": "", "journal": ""}]

GALLERY IMAGES

