


GOAT ANTI-CFD / ADIPSIN (AA87-96) ANTIBODY

SKU: EB12706



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	CFD complement factor D (adipsin) ADIPSIN ADN DF PFD C3 convertase activator D component of complement (adipsin) complement factor D complement factor D preproprotein properdin factor D
Accession ID	NP_001919.2
Blocking Peptide	EBP12706
Immunogen	Peptide with sequence C-QPEPSKRLYD, from the internal region of the protein sequence according to NP_001919.2.
Peptide Sequence	C-QPEPSKRLYD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Pig, Cow
Reactive Species	Human
Human Gene ID	1675
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 30kDa band observed in Human Adipose lysates (calculated MW of 27kDa according to NP_001919.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

