

GOAT ANTI-CHD1 (C TERMINUS) ANTIBODY

SKU: EB07840



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	DKFZp686E2337 chromodomain helicase DNA binding protein 1 CHD1
Accession ID	NP_001261.2
Blocking Peptide	EBP07840
Immunogen	Peptide with sequence C-EHKSTPEHTWSSRK, from the C Terminus of the protein sequence according to NP_001261.2.
Peptide Sequence	C-EHKSTPEHTWSSRK
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	1105
Mouse Gene ID	12648
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments in Human Liver, Lung and Skeletal Muscle lysates gave no specific signal but low background (at antibody concentration up to 0.3µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA