



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CHN1 ANTIBODY

SKU: EB05206

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Chimerin 1 (GTPase-activating protein, rho, 2)|n-chimerin|RhoGAP2|RHOGAP2|ARHGAP2|CHN|chimerin (chimaerin) Alias

Names

1|CHN1

Accession

NP 001813

Immunogen Peptide with sequence C-ELLIKNEDILF, from the C Terminus of the protein sequence according to NP_001813.

Peptide

C-ELLIKNEDILF

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Instructions Refrigerated

Predicted Species

Human, Rat, Monkey

Human

1123 **Gene ID**

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Western

Western blot: No signal obtained yet but low background observed in Mouse Brain extracts at upto 2µg/ml. Blot