

GOAT ANTI-CHN1 ANTIBODY

SKU: EB05206



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Chimerin 1 (GTPase-activating protein, rho, 2) n-chimerin RhoGAP2 RHOGAP2 ARHGAP2 CHN chimerin (chimaerin) 1 CHN1
Accession ID	NP_001813
Immunogen	Peptide with sequence C-ELLIKNEDILF, from the C Terminus of the protein sequence according to NP_001813.
Peptide Sequence	C-ELLIKNEDILF
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Rat, Monkey
Human Gene ID	1123
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Mouse Brain extracts at upto 2µg/ml.