






GOAT ANTI-CHRNB2 ANTIBODY

SKU: EB08480



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Instructions

Synonym / Alias

neuronal nicotinic acetylcholine receptor beta 2|cholinergic receptor, nicotinic, beta polypeptide 2 (neuronal)|nAChRB2|EFNL3|cholinergic receptor, nicotinic, beta 2 (neuronal)|CHRNB2

Names Accession

NP_000739.1

Blocking

ID

EBP08480

Peptide

Immunogen

Peptide with sequence QPRHHCARQRLR, from the internal region of the protein sequence according to NP_000739.1.

Peptide Sequence

QPRHHCARQRLR

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human, Mouse, Rat, Dog, Pig, Cow

Species Reactive

Human

Species Human

Gene ID

1141

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results Paraffin embedded Human Uterus. Recommended concentration: 2.5µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:64000.

Limit

Western Appr Blot Appr

Approx. 40-48 kDa band observed in Human Brain (Cerebellum),) lysates (calculated MW of 57.0 kDa

according to NP_000739.1). This molecular weight is routinely observed by other sources. Recommended $\frac{1}{2}$

concentration: 1-3µg/ml. Primary incubation 1 hour at room temperature.

Application

Type

Pep-ELISA, WB, IHC

GALLERY IMAGES



Telephone: (650) 697-3600





