

GOAT ANTI-CLCA1 (AA872-884) ANTIBODY

SKU: EB08583



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias CaCC|chloride channel, calcium activated, family member 1|CLCA1

Names

Accession ID NP_001276.2

Blocking Peptide EBP08583

Immunogen Peptide with sequence PETPSPDETSAPC, from the internal region of the protein sequence according to NP_001276.2.

Peptide Sequence PETPSPDETSAPC

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 1179

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Limit Antibody detection limit dilution 1:64000.

Western Blot Approx. 100kDa band observed in Human Kidney and Lung lysates (calculated MW of 100kDa according to NP_001276.2). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

