

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-CLCA1 (AA872-884) ANTIBODY

SKU: EB08583



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

CaCC|chloride channel, calcium activated, family member 1|CLCA1

Names

ID

Accession

NP_001276.2

Blocking

EBP08583 **Peptide**

Peptide with sequence PETPSPDETSAPC, from the internal region of the protein sequence according to **Immunogen**

NP 001276.2.

Peptide

PETPSPDETSAPC Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

1179 **Gene ID**

Product

Grade **ELISA**

Detection Antibody detection limit dilution 1:64000.

Limit

Western Approx. 100kDa band observed in Human Kidney and Lung lysates (calculated MW of 100kDa according to

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Blot NP_001276.2). Recommended concentration: 0.3-1μg/ml. Primary incubation was 1 hour.

Application

Pep-ELISA, WB Type

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa