





## **GOAT ANTI-CLEC16A (AA1040-53) ANTIBODY**

**SKU:** EB09210



Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

CLEC16A|C-type lectin domain family 16, member A|Gop-1|KIAA0350|MGC111457

Alias Names

ID

**Accession** 

NP\_056041.1

**Blocking** 

EBP09210 **Peptide** 

Peptide with sequence EEAACAEPVGTAED, from the C Terminus of the protein sequence according to **Immunogen** 

NP 056041.1.

**Peptide** 

Sequence

**EEAACAEPVGTAED** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive

Human **Species** 

Human

23274 Gene ID

**Product** Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:1000.

Limit

Western Approx. 110kDa band observed in Human Kidney and Liver lysates (calculated MW of 118kDa according to

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**Blot** NP\_056041.1). Recommended concentration: 1-3µg/ml.

**Application** 

Pep-ELISA, WB Type

## **GALLERY IMAGES**



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa