

GOAT ANTI-CLEC16A (AA1040-53) ANTIBODY

SKU: EB09210



250kDa
150kDa
100kDa
75kDa

50kDa

37kDa

25kDa
20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	CLEC16A C-type lectin domain family 16, member A Gop-1 KIAA0350 MGC111457
Accession ID	NP_056041.1
Blocking Peptide	EBP09210
Immunogen	Peptide with sequence EEAACAEPVGTAED, from the C Terminus of the protein sequence according to NP_056041.1.
Peptide Sequence	EEAACAEPVGTAED
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	23274
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx. 110kDa band observed in Human Kidney and Liver lysates (calculated MW of 118kDa according to NP_056041.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

