


GOAT ANTI-CLO29 (AA30-43) ANTIBODY

SKU: EB11549



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	C12orf29 chromosome 12 open reading frame 29 CLO29 DKFZp313K0436 DKFZp434N2030 DKFZp686L04169 FLJ38158 hypothetical protein LOC91298 MGC102978
Accession ID	NP_001009894.2
Blocking Peptide	EBP11549
Immunogen	Peptide with sequence C-QPFKVLATETVSHK, from the internal region of the protein sequence according to NP_001009894.2.
Product Comments	This antibody is expected to recognize isoform 1 (Q8N999-1) and isoform 3 (Q8N999-3) as described in UniprotKB.
Peptide Sequence	C-QPFKVLATETVSHK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Pig, Sheep, Cow
Reactive Species	Human
Human Gene ID	91298
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Approx 38kDa band observed in lysates of cell line HeLa (calculated MW of 37.5kDa according to NP_001009894.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

