

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600


GOAT ANTI-CLUSTERIN / APOJ (MOUSE) ANTIBODY


SKU: EB08920



Telephone: (650) 697-3600









SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

testosterone repressed prostate message|complement lysis inhibitor|Apolipoprotein Alias

J|Sugp-2|Sgp2|Sgp-2|SP-40|D14Ucla3|Cli|Al893575|Clu|clusterin|Apo| Names

Usage Immunofluorescence: Strong expression of the protein seen in the membranes and

cytoplasm of PC12 cells. Recommended concentration: 10µg/ml. Summary

Accession

NP_038520.2

Blocking

ID

EBP08920 **Peptide**

Peptide with sequence C-EKALQEYRRKSRAE, from the C Terminus of the protein sequence according to **Immunogen**

NP 038520.2.

Peptide

Sequence

C-EKALQEYRRKSRAE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Mouse, Rat

Reactive

Mouse, Rat **Species**

Mouse

Gene ID

12759

Rat Gene ID 24854

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:64000.

Limit

Approx. 38kDa band observed in Mouse and Rat Brain lysates.(calculated MW of 51.7kDa according to Mouse NP_038520.2 and 51.4kDa according to Rat XP_006252156.1). This molecular weight is routinely observed by

Western **Blot**

other sources and thought to correspond to the alpha chain of APOJ. Recommended concentration:

0.1- $0.3\mu g/ml$. Primary incubation 1 hour at room temperature.

Application

Type

Pep-ELISA, WB, IF







GALLERY IMAGES

