

GOAT ANTI-COILIN (AA73-86) ANTIBODY

SKU: EB11577

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias CLN80|COIL|coilin|coilin p80|p80-coilin

Names

Accession ID NP_004636.1

Blocking Peptide EBP11577

Immunogen Peptide with sequence RLVRDNDCLRVKLE, from the internal region of the protein sequence according to NP_004636.1.

Peptide Sequence RLVRDNDCLRVKLE

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Rat, Cow

Reactive Species Human

Human Gene ID 8161

Rat Gene ID 50998

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Limit Antibody detection limit dilution 1:128000.

Western Blot Approx 80kDa band observed in nuclear lysates of cell line MCF7 (calculated MW of 62.6kDa according to NP_004636.1). Recommended concentration: 1-3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

