


GOAT ANTI-COL4A3 (AA341-353) ANTIBODY

SKU: EB10708



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Instructions	
Synonym / Alias Names	tumstatin collagen alpha-3(IV) chain collagen IV, alpha-3 polypeptide collagen, type IV, alpha 3 (Goodpasture antigen) COL4A3
Accession ID	NP_000082.2
Blocking Peptide	EBP10708
Immunogen	Peptide with sequence C-RGPTEYYDTYQEK, from the internal region of the protein sequence according to NP_000082.2.
Peptide Sequence	C-RGPTEYYDTYQEK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated

Predicted Species	Human, Dog, Cow
Reactive Species	Human
Human Gene ID	1285
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Kidney shows staining at the luminal side of the epithelial cell membranes. Recommended concentration, 4-6µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

