


GOAT ANTI-COL4A3BP (AA396-411) ANTIBODY

SKU: EB08911



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	lipid-transfer protein CERTL ceramide transporter alpha 3 type IV collagen binding protein StAR-related lipid transfer (START) domain containing 11 START domain containing 11 Goodpasture antigen-binding protein STARD11 GPBP FLJ20597 CERTL CERT collagen, type IV, alpha 3 (Goodpasture antigen) binding protein COL4A3BP
Accession ID	NP_005704.1; NP_112729.1
Blocking Peptide	EBP08911
Immunogen	Peptide with sequence C-KREDSWQKRLDKETEK, from the internal region of the protein sequence according to NP_005704.1; NP_112729.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_005704.1; NP_112729.1).
Peptide Sequence	C-KREDSWQKRLDKETEK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Cow
Reactive Species	Human
Human Gene ID	10087
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx. 70kDa band observed in Human Heart lysates (calculated MW of 70.8kDa according to NP_005704.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

